

HAZMAX Hazardous Containment Work Area

HazMax Floor-Mounted Hoods are designed to safely isolate hazardous materials and procedures from the work area. Hood sizes include 76" and 100" wide models in 48" and 72" depths. Ideal for enclosing potentially hazardous mixing and dispensing procedures or for providing a convenient safe area to store hazardous materials.



HAZMAX Floor Mount Hood is shown with optional service features and accessories.

HAZMAX Floor-Mounted Hood Features

Hood features a chemical resistant, non-sparking construction. Upper and lower air to bypass inlets together with an advanced baffle design maintain a uniform airflow through the fume chamber. Horizontal slide doors are of aluminum framed tempered glass. Optional door systems are available.

Electrical Components

Hood is equipped with vapor proof or explosion proof lighting. All electrical components and services are UL and CSA listed. Other voltage Hertz and phase options are available for international and industry requirements.

Corrosion Resistant Construction

The molded one-piece containment basin safely holds any accidental spills of hazardous liquids or materials. The chemical resistant high strength, raised fiberglass grating supports equipment and personnel. Non-skid grate decking is removable for cleaning and decontamination. Outlet connection is available for drainage basin.

Note: There may be procedures that require personnel entry into floor-mounted hood. Person should have protection/respirators. Consult laboratory safety officer for correct personnel protection equipment.

