



CAT. NO. 21421
(WITH OPTIONAL BASE CABINET
AND WORKSURFACE)

CAT. NO. 21521
(WITH OPTIONAL TABLE
AND WORKSURFACE)

CAT. NO. 21621
(WITH OPTIONAL ACID/BASE CABINETS
AND WORKSURFACE)

CAT. NO. 21821
(WITH OPTIONAL HORIZONTAL
SLIDING SASH, BASE CABINETS
AND WORKSURFACE)

SECTION A-A

OPTIONAL EQUIPMENT AND SERVICES

1. PLUMBING SERVICE
2. ELECTRICAL SERVICES
3. EXHAUST SYSTEMS
4. FIRE EXTINGUISHER
5. AIR FLOW ALARM
6. HORIZONTAL SLIDING SASH
7. COMBINATION SASH

NOTE:
ALL ELECTRICAL COMPONENTS
UL AND CSA LISTED.

AUXAIRFAM21X21-SV

UNIFLOW AUXILIARY AIR FUME HOOD			
SCALE: FULL	APPROVED BY:	DATE:	REVISION:
DRAWN BY: TLF	REP:	QUOTE No:	DATE:
DATE: 08/22/03	PROJECT MGR:	DATE:	DATE:
HEMCO Corporation		© PLANS AND DRAWINGS ARE COPYRIGHTED	DWG. NO. M01-3270
111 POWELL RD. INDEPENDENCE, MO. 64056		PH:816-796-2900 FAX:816-796-3333	DATE: _____ SHEET. 1 OF 1