



CAT. NO. 17041
(WITH OPTIONAL BASE CABINET
AND WORKSURFACE)

CAT. NO. 17051
(WITH OPTIONAL TABLE,
WORKSURFACE, AND
FULL SERVICES.)

CAT. NO. 17061
(WITH OPTIONAL ACID/BASE CABINETS
AND WORKSURFACE)

CAT. NO. 17081
(WITH OPTIONAL BASE CABINETS
HORIZONTAL SLIDING SASH
AND WORKSURFACE)

SECTION A-A

OPTIONAL EQUIPMENT AND SERVICES

1. PLUMBING SERVICE
2. ELECTRICAL SERVICES
3. EXHAUST SYSTEMS
4. STATIC PRESSURE MONITOR
5. AIR FLOW ALARM
6. HORIZONTAL SLIDING SASH
7. COMBINATION SASH

NOTE:
ALL ELECTRICAL COMPONENTS
UL AND CSA LISTED.

CLNAIRFAM17XX1-SV

UNIFLOW CLEAN AIRE I FUME HOOD

SCALE: 1"=24"	APPROVED BY:	DATE:	REVISION:	DATE:
DRAWN BY: TLF	REP:	QUOTE No:	XXXX XXXX	XX/XX/XX
DATE: 08/25/03	PROJECT MGR:	DATE:		
HEMCO Corporation www.HEMCOcorp.com 111 POWELL RD INDEPENDENCE, MO 64056 PH:816-796-2900 FAX:816-796-3333			© PLANS AND DRAWINGS ARE COPYRIGHTED	DWG. NO. M01-3275
			BY: _____	SHEET. 1 OF 1
			DATE: _____	