



**MODEL 21401**  
(WITH OPTIONAL BASE CABINET AND WORKSURFACE)

**MODEL 21501**  
(WITH OPTIONAL TABLE AND WORKSURFACE)

**MODEL 21601**  
(WITH OPTIONAL ACID/BASE CABINET AND WORKSURFACE)

**MODEL 21801**  
(WITH OPTIONAL HORIZONTAL SLIDING SASH, BASE CABINETS AND WORKSURFACE)

**SECTION A-A**

**OPTIONAL EQUIPMENT AND SERVICES**

1. PLUMBING SERVICE
2. ELECTRICAL SERVICES
3. EXHAUST SYSTEMS
4. FIRE EXTINGUISHER
5. AIR FLOW ALARM
6. HORIZONTAL SLIDING SASH
7. COMBINATION SASH

**NOTE:**  
ALL ELECTRICAL COMPONENTS  
UL AND CSA LISTED.

CONVENTFAM21X01-SV

UNIFLOW CONVENTIONAL/VAV FUME HOODS				
SCALE: FULL	APPROVED BY:	DATE:	REVISION:	DATE:
DRAWN BY: DON	REP:	QUOTE No:	XXXX XXXX	XX/XX/XX
DATE: 09/27/04	PROJECT MGR:	DATE:		
<b>HEMCO Corporation</b>		© PLANS AND DRAWINGS ARE COPYRIGHTED		DWG. NO. M01-3184
111 POWELL RD INDEPENDENCE, MO. 64056 PH:816-796-2900 FAX:816-796-3333		BY _____ DATE _____		SHEET. 1 OF 1