



CAT. NO. 13092
(WITH OPTIONAL BASE CABINET)
AND WORKSURFACE

CAT. NO. 13692
(WITH OPTIONAL TABLE
AND WORKSURFACE)

CAT. NO. 14892
(WITH OPTIONAL ACID CABINET
AND WORKSURFACE)

OPTIONAL EQUIPMENT AND SERVICES

1. PLUMBING SERVICE
2. ELECTRICAL SERVICES
3. EXHAUST SYSTEMS
4. FIRE EXTINGUISHER
5. AIR FLOW ALARM
6. HORIZONTAL SLIDING SASH

NOTE:
ALL ELECTRICAL COMPONENTS
ARE UL & CSA LISTED.

UNIFLOW CE LABORATORY FUME HOOD SERIES 02 30", 36" & 48" WIDTHS			
SCALE:	APPROVED BY:	DATE:	REVISION:
DRAWN BY: DON	REP:	QUOTE No:	B- UPDATED TO CURRENT SPEC. 09/18/06
DATE: 08/19/03	PROJECT MGR:	DATE:	A- UPDATED TO CURRENT SPEC. 05/03/05
HEMCO Corporation www.HEMCOcorp.com 111 POWELL RD INDEPENDENCE, MO 64050 PH:816-796-2900 FAX:816-796-3333			DWG. NO. M01-3255
BY R. E. Hill DATE 08/19/03			SHEET. 1 OF 1