



CAT. NO. 13411  
(WITH OPTIONAL BASE CABINET)

CAT. NO. 13511  
(WITH OPTIONAL TABLE)

CAT. NO. 13611  
(WITH OPTIONAL ACID/BASE CABINETS)

CAT. NO. 13811  
(WITH OPTIONAL BASE CABINET)

SECTION A-A

OPTIONAL EQUIPMENT AND SERVICES

1. PLUMBING SERVICE
2. ELECTRICAL SERVICES
3. EXHAUST SYSTEMS
4. AIR FLOW ALARM

**NOTE:**  
THE FUME CHAMBER, VARIFLOW BAFFLE SYSTEM, INTEGRAL DISHED WORKSURFACE, AND THE DRAINAGE TROUGH ARE FABRICATED FROM PVC.

**NOTE:**  
ALL ELECTRICAL COMPONENTS ARE UL & CSA LISTED.

ACIDDIGFAM13X11-SV

### UNIFLOW ACID DIGEST FUME HOOD SERIES 01 WITH VAPOR PROOF LIGHT

SCALE: 1" = 16"	APPROVED BY:	DATE:	REVISION:	DATE:
DRAWN BY: DON	REP:	QUOTE No:		
DATE: 12/02/06	PROJECT MGR:	DATE:		
<b>HEMCO Corporation</b>		© PLANS AND DRAWINGS ARE COPYRIGHTED		DWG. NO. M01-4160
www.HEMCOcorp.com 111 POWELL RD INDEPENDENCE, MO 64056 PH:816-796-2900		BY <i>R. E. Hill</i> DATE 12/02/06		SHEET. 1 OF 1