



CAT. NO. 45411
(WITH OPTIONAL BASE CABINET AND WORKSURFACE)

CAT. NO. 45511
(WITH OPTIONAL TABLE AND WORKSURFACE)

CAT. NO. 45611
(WITH OPTIONAL ACID/BASE CABINETS AND WORKSURFACE)

CAT. NO. 45811
(WITH OPTIONAL HORIZONTAL SLIDING SASH AND BASE CABINETS AND WORKSURFACE)

SECTION A-A

OPTIONAL EQUIPMENT AND SERVICES

1. PLUMBING SERVICE
2. ELECTRICAL SERVICES
3. EXHAUST SYSTEMS
4. FIRE EXTINGUISHER
5. AIR FLOW ALARM
6. HORIZONTAL SLIDING SASH
7. COMBINATION SASH

NOTE:
ALL ELECTRICAL COMPONENTS
UL AND CSA LISTED

UNIFLOW SE LABORATORY FUME HOOD SERIES 45 WITH VAPOR PROOF LIGHT			
SCALE: 1" = 24"	APPROVED BY:	DATE:	REVISION:
DRAWN BY: DON	REP:	QUOTE No:	DATE:
DATE: 8/22/06	PROJECT MGR:	DATE:	
HEMCO Corporation www.HEMCOcorp.com 111 POWELL RD INDEPENDENCE, MO 64056 PH:816-796-2900 FAX:816-796-3333		PLANS AND DRAWINGS ARE COPYRIGHTED BY <i>R. E. Hill</i> DATE 08/22/06	DWG. NO. M01-4099 SHEET. 1 OF 1

BYPASSFAM 45X11-SV