



**NOTE:**  
CONSTRUCTION IS OF NON SPARKING MATERIALS THROUGHOUT.

	A	B	C	BASIN CAPACITY	CATALOG NUMBER
1	76"	48"	40"	42 GAL	50864
2	76"	72"	64"	70 GAL	50866
3	100"	48"	40"	59 GAL	50844
4	100"	72"	64"	88 GAL	50886

HAZMAXCAT-SV

HAZMAX HAZARDOUS CONTAINMENT HOOD VAPOR PROOF LIGHTING				
SCALE: 1"=20"	APPROVED BY:	DATE:	REVISION:	DATE:
DRAWN BY: DON	REP:	QUOTE No:	XXXX XXXX	XX/XX/XX
DATE: 09/20/04	PROJECT MGR:	DATE:		
<b>HEMCO Corporation</b>			© PLANS AND DRAWINGS ARE COPYRIGHTED	DWG. NO. H19-3664
111 POWELL RD INDEPENDENCE MO. 64056 PH:816-796-2900 FAX:816-796-3333			BY: _____ DATE: _____	SHEET. 1 OF 1