



MODEL 12401
(WITH OPTIONAL BASE CABINET)

MODEL 12501
(WITH OPTIONAL TABLE)

MODEL 12601
(WITH OPTIONAL ACID/BASE CABINETS)

SECTION A-A

OPTIONAL EQUIPMENT AND SERVICES

1. PLUMBING SERVICE
2. ELECTRICAL SERVICES
3. EXHAUST SYSTEMS
4. AIR FLOW ALARM

NOTE:
THE FUME CHAMBER, VARIFLOW BAFFLE SYSTEM, INTEGRAL DISHED WORKSURFACE, AND THE DRAINAGE TROUGH ARE FABRICATED FROM 316 SERIES STAINLESS STEEL.

NOTE:
ALL ELECTRICAL COMPONENTS ARE UL & CSA LISTED.

PERCACIDFAM12X01-SV

UNIFLOW PERCHLORIC ACID FUME HOOD SERIES 01 VAPOR PROOF LIGHT & WASH DOWN SYSTEM			
SCALE: FULL	APPROVED BY:	DATE:	REVISION:
DRAWN BY: DON	REP:	QUOTE No:	XXXX XXXX
DATE: 09/08/03	PROJECT MGR:	DATE:	XX/XX/XX
HEMCO Corporation		© PLANS AND DRAWINGS ARE COPYRIGHTED	DWG. NO. M01-3282
111 POWELL RD. INDEPENDENCE, MO. 64056 PH:816-796-2900	FAX:816-796-3333	BY: _____ DATE: _____	SHEET. 1 OF 1