

Fume Hood Work Surfaces

Fume Hood Work Surfaces are available to suit the specific needs of your Lab, ranging from work surfaces built to withstand heavy acid solvent use, to light weight general chemical resistant use. All HEMCO worksurfaces are tested to U.L. performance / high heat requirement.




U.L 1805 Certified

Recommended for general laboratory work where high concentrations of solvent, acids, and alkalis, are in use.

Solid Epoxy Resin Fume Hood Work Surfaces

Epoxy Resin Work surfaces are recommended for general laboratory work where high concentrations of solvent, acids, and alkalis are being worked with. Work surface is made of cast thermosetting epoxy resin, creating a molded solid 1-1/4" thick top which is extremely hard, and resistant to abrasions and heat. Surface is dished 3/8" to contain spillage. Designed to fit HEMCO Fume Hoods. Color: Black

Size	Cat. No.	Size	Cat. No.
36"W x 30"D	20315	36"W x 36"D	20325
48"W x 30"D	20415	48"W x 36"D	20425
60"W x 30"D	20515	60"W x 36"D	20525
72"W x 30"D	20615	72"W x 36"D	20625
96"W x 30"D	20815	96"W x 36"D	20825



U.L 1805 Certified

Recommended for general Laboratory work where greater chemical resistance is required.

Solid Phenolic Resin Fume Hood Work Surfaces

Phenolic resin work surfaces are recommended for general laboratory work where greater chemical resistance is required. Work surface combines excellent chemical resistance, with superior scratch and wear resistance. Work surface is 1-1/4" thick and is dished 1/4". Color: Black

Size	Cat. No.	Size	Cat. No.
36"W x 30"D	21318	36"W x 36"D	22318
48"W x 30"D	21418	48"W x 36"D	22418
60"W x 30"D	21518	60"W x 36"D	22518
72"W x 30"D	21618	72"W x 36"D	22618
96"W x 30"D	21818	96"W x 36"D	22818



U.L 1805 Certified

Recommended for general Laboratory work where chemical resistance is required.

Solid Composite Resin Fume Hood Work Surfaces

Dura-Lite composite resin work surfaces are recommended for general laboratory work where chemical resistance is required. Work surface is reinforced with advanced composite materials 1-1/4" thick, work area is dished 3/8" to contain spillage. Designed to fit HEMCO Fume Hoods. Color: Gray

Size	Cat. No.
36"W x 30"D	20313
48"W x 30"D	20413
60"W x 30"D	20513
72"W x 30"D	20613
96"W x 30"D	20813

Visit www.HEMCOcorp.com

