
## Canopy Hoods in Fiberglass or Stainless Steel

Canopy hoods are designed to collect and exhaust corrosive vapors, heat, steam, and odors when mounted over areas with water baths, hot plates, or portable equipment.

Fiberglass Wall Canopy



Fiberglass Canopy hoods are designed to collect and exhaust corrosive vapors, heat, steam, and odors when mounted over areas with water baths, hot plates, or portable equipment. Standard color: Silver Beige



Side Wall Panels (fiberglass or stainless steel) kit includes one wall panel and hardware Cat. No. 23000

Vapor Proof Light


Rear Wall Panel (consult factory for details)

Blower Switch 115V/15A

& Switch Cat. No. 50035	Cat. No. 50026				
Width A	36"	48"	72"	96"	Side Wall Panels
Fiberglass Cat. No.	13030	13040	13060	13080	23000
Fiberglass Rear Wall Panel	23016	23026	23036	23046	
Diameter D	6"	8"	10"	(2) 10"	

Fiberglass Island Canopy






Suspended Island Canopy Hood Mounting Kit Includes all components to secure to ceiling including 6' long threaded rods. Cat. No. 23001

Side Wall Panels (fiberglass or stainless steel) kit includes one wall panel and hardware Cat. No. 23000 **Rear Wall Panel** 


Vapor Proof Light & Switch Cat. No. 50035 (consult chart for Cat.No.)

	(contain chart for Caurto)				
Width A	36"	48"	72"	96"	Side Wall Panels
Fiberglass Cat. No.	23030	23040	23060	23080	23000
Fiberglass Rear Wall Panel	23016	23026	23036	23046	-
Diameter D	6"	8"	10"	(2) 10"	

Stainless Steel Wall Canopy



Stainless Steel Canopy hoods are constructed of type 304 stainless steel and designed to exhaust heat, steam when mounted over water baths, hot plates, or portable equipment. In sizes 36", 48", 72", 96" wide.



Side Wall Panels (fiberglass or stainless steel) kit includes one wall panel and hardware Cat. No. 23011


Vapor Proof Light & Switch Cat. No. 50035

**Rear Wall Panel** (consult factory for details) Blower Switch 115V/15A Cat. No. 50028

Width A	36"	48"	72"	96"	Side Wall Panels
Stainless Steel Cat No.	13130	13140	13160	13180	23011
Stainless Steel Rear Wall Panel	23116	23126	23136	23146	
Diameter D	6"	8"	10"	(2) 10"	

Stainless Steel Island Canopy





Suspended Island Canopy Hood Mounting Kit Includes all components to secure to ceiling including 6' long threaded rods. Cat. No. 23001

Side Wall Panels (fiberglass or stainless steel) kit includes one wall panel and hardware Cat. No. 23011


Blower Switch 115V/15A (blower not included) Cat. No. 50027-2

1101 111014404 / 0441 1101 00021 2				· -	
Width A	36"	48"	72"	96"	Side Wall Panels
Stainless Steel Cat No.	23130	23140	23160	23180	23011
Stainless Steel Rear Wall Panel	23116	23126	23136	23146	
Diameter D	6"	8"	10"	(2) 10"	

# Canopy Hoods in Fiberglass or Stainless Steel Canopy hoods an be used for storage of chemicals, freeing up space in fume hoods.







Corner Canopy Hood & Wall Mounting Kit Includes all components to secure the canopy hood to a wall. Cat. No. 13000

Corner Composite Resin Worksurface Comes dished to contain spillage. Cat. No. 13503

Corner Table Features welded steel construction Cat. No. 13502

Rear Wall Panel (2) required Qty (1) 36" high X 45" wide Cat. No. 23056

Side Wall Panels (fiberglass) 36" h X 25" w, kit includes one wall panel and hardware Cat. No. 23006

Blower Switch 115V/15A Cat. No. 50028

#### **Custom Fiberglass Canopy Hoods**



Canopy with (2) side walls, for access from both sides.



Canopy with wall mounting brackets.



Corner Canopy with light, side walls and corner table unit.




Island Canopy with ducting and work table.

#### Corner Canopy (Stainless Steel)







Corner Canopy Hood & Wall Mounting Kit Includes all components to secure the canopy hood to a wall. Cat. No. 13000

Corner Composite Resin Worksurface Comes dished to contain spillage. Cat. No. 13503

Corner Table Features welded steel construction Cat. No. 13502

### **Custom Stainless Steel Canopy Hoods**



Canopy with (2) side walls, for access from both sides



Side Wall Panels (stainless steel)

36" h X 25" w, kit includes one wall

Blower Switch 115V/15A


Cat. No. 50027-2

panel and hardware Cat. No. 23004

Canopy with rear and side walls, worksurface, and adjustable table.



Canopy with Plexiglas sides, and hinged



48" wall canopy, with light switch.







