## **VSE Balance Enclosure**

Versatile vented balance enclosure is ideal for critical weighing procedures involving hazardous powders and liquids.



VSE - Cat. No. 84800 shown above with optional Mobile Table and Phenolic Worksurface.

DESCRIPTION CATALOG NUMBER			
	24" W	36"W	48"W
VSE Balance Enclosure 24" Depth X 22" height	82400	83600	84800
Phenolic Work surface	12418	13618	14818
Mobile Table		50318	50418
Flex Ducting 6" Dia. X 12' long	80052	80052	80052
Exhaust Fan for remote mount 115 volts, 60Hz	51051	51052	51053

- VSE comes in 24", 36" and 48" widths to accommodate all analytical balances.
- Efficient air flow design with air foils and by pass directs airborne particulate to baffled exhaust thereby providing superior air flow and containment performance for user protection.
- Ergonomic design provides added user comfort.
- Constructed of chemical resistant metal framing with 1/4" thick clear tempered safety glass side panels and viewing sash.
- The viewing sash is angled 15 degrees for ease of viewing comfort with 8" reach in opening height.
  sash swings up to provide 20" of access opening.
- 2 service ports are located in lower right and left rear wall.



VSE - Cat. No. 83600 shown above with optional Mobile Table & Phenolic Worksurface.





