


Ductless Demonstration Hood

The HEMCO Ductless Demo Hood effectively removes fumes and odors while providing a clear, unobstructed view of the interior from all sides. The demonstration hood is designed for classroom and group observation situations.



- A Class DH-1 ductless hood workstation Equipped with filtration devices designed to control non-toxic chemicals, nuisance odors and particulate.
- Dimensions 35" wide X 24" deep X 30.5 " high with phenolic resin worksurface attached.
- Viewing height is 24", through hinged clear acrylic viewing sash, shatterproof side and rear panels.
- Equipped with a vapor proof T-5 fluorescent light fixture, and control switch.
- Front panel includes light switch, 158 CFM variable speed blower, switch and filter service light, 8' Power cord.
- Activated Carbon Filtration for particulate, fumes, odors, non-hazardous chemical vapors.



Specifications:

Demo Hood: Cat. No. 53015 115 V AC, 60 Hz
International Electrical Cat. No. 53015 230 V AC, 50 Hz

Mobile Table: Cat. No. 50318 37" X 25"
Height adjustable from 30" to 36"

Non-Adjustable Table Cat. No. 50315
Fixed Height 36"

Replacement Filter: Cat. No. 52005 Qty 2 pack

Shipping Weight: 250 lbs. or 113.64 kilos.

Visit www.HEMCOcorp.com



HEMCO Corporation
711 S. Powell Road
Independence, MO 64056

HEMCO

Phone (816) 796-2900
Fax (816) 796-3333
DDH-08.11