

Emergency Shower Decontamination Booth

The HEMCO Emergency Shower provides a complete solution to your emergency shower and decontamination needs. Providing a immediate full bodied shower & eye/face wash.

- ✓ Full body shower and eye/face wash, both with immediate response activation handle.
- ✓ Fully Factory Plumbed to a single inlet, 1" diameter pipe. 30 GPM water supply required.
- ✓ Durable one piece fiberglass construction with glass smooth white surfaces and covered corners for ease of clean-up. Non-skid surface with self draining containment basin, and 3" diameter drain outlet.
- ✓ Emergency Shower Decontamination Booth is pre plumbed to single point 1/4 "IPT. connection. Meets ANSI and OSHA requirements. Tempered warm water supply connection, is recommended.

Specifications:

Emergency Drench Shower Booth- molded of white fiberglass and features seamless construction. Unit is equipped with pull rod activated overhead drench shower and rear wall mounted eye/face wash **Cat No. 16601**



Accessories:

1. **Translucent / Frosted Strip Curtain** - provides for easy access into booth, while containing contaminated water splash. **Cat No. 16650**
2. **Alarm Horn**- audible warning sound when system is actuated. **Cat No. 16680**
3. **Alarm Light**- Visual warning light when system is activated. **Cat No. 16670**
4. **Hand Held Body wash**- Squeeze valve nozzle on 8' coiled hose. Can be used to rinse all body parts while victim is sitting or lying down. **Cat No. 16005**
5. **Grab Bars**- Support victim and rescue personnel during decontamination. **Cat No. 16620**



Visit www.HEMCOcorp.com