

Fume Hood Fire Extinguishers

HEMCO Fume hood Fire Extinguishers, for use in situations where there is a possibility of exposure to flammables. In addition to the peace of mind, Extinguishers include a three-year limited warranty.

NFPA 45, paragraph 6.10.2 in accordance with the Standard on Fire Protection for Laboratories Using Chemicals, states that spaces such as Fume Hood Enclosures be equipped with fire extinguishers. HEMCO offers 2.7, 5.4 and 10lb Clean Agent FE-241 fire extinguishers. A fuseable link mounted inside the fume chamber senses excessive heat and automatically discharges the system. After activation the fire is extinguished keeping both users and the Lab environment safe from any further threat of fire. These units are factory installed when ordered with a hood, or they can be field mounted on most pre-existing hoods.



Remote System Status Light



Remote Discharge Alarm



Front Mount Manual Discharge Cable



Specifications / Catalog Numbers

- 2.7 lb. Extinguisher** **Cat No. 53050**
 Easily mounts in existing fume hoods measuring up to 4'
- 5.4 lb. Extinguisher** **Cat No. 53125**
 Easily mounts in existing fume hoods measuring up to 8'
- 10 lb. Extinguisher** **Cat No. 53200**
 Easily mounts in existing fume hoods measuring over 8'
- Remote System Status Light** **Cat No. 53006**
 The green warning light provides constant monitoring of the system
- Remote (audible) Discharge Alarm** **Cat No. 53007**
 Features integrated alarm that provides around the clock discharge alert
- Front Mount Discharge Cable** **Cat No. 53008**
 Discharge cable kit includes cable, escutcheon plate, pull handle, safety pin and safety tie

Visit www.HEMCOcorp.com