

# Fume Hood Work Surfaces

Fume Hood Work Surfaces are available to suit the specific needs of your Lab, ranging from work surfaces built to withstand heavy acid solvent use, to light weight general chemical resistant use. All HEMCO worksurfaces are tested to U.L. performance / high heat requirement.



**Recommended for general laboratory work where high concentrations of solvent, acids, and alkalis, are in use.**

## Solid Epoxy Resin Fume Hood Work Surfaces



U.L 1805 Certified

Epoxy Resin Work surfaces are recommended for general laboratory work where high concentrations of solvent, acids, and alkalis are being worked with. Work surface is made of cast thermosetting epoxy resin, creating a molded solid 1 1/4" thick top which is extremely hard, and resistant to abrasions and heat. Surface is dished to contain spillage. Designed to fit HEMCO Fume Hoods. Color: Black

Epoxy Resin Work Surface Sizes		
Wide	Deep	Cat. No.
36"	30"	20315
48"	30"	20415
60"	30"	20515
72"	30"	20615
96"	30"	20815



**Recommended for general Laboratory work where greater chemical resistance is required.**

## Solid Phenolic Resin Fume Hood Work Surfaces



U.L 1805 Certified

Phenolic resin work surfaces are recommended for general laboratory work where greater chemical resistance is required. Work surface combines excellent chemical resistance, with superior scratch and wear resistance. Work surface is 1" thick and is flat. Designed to fit HEMCO Fume Hoods. Color: Black

Phenolic Resin Work Surface Sizes		
Wide	Deep	Cat. No.
36"	30"	20318
48"	30"	20418
60"	30"	20518
72"	30"	20618
96"	30"	20818



**Recommended for general Laboratory work where chemical resistance is required.**

## Solid Composite Resin Fume Hood Work Surfaces



U.L 1805 Certified

Dura-Lite composite resin work surfaces are recommended for general laboratory work where chemical resistance is required. Work surface is reinforced with advanced composite materials 1 1/4" thick, work area is dished 3/8" to contain spillage. Designed to fit HEMCO Fume Hoods. Color: Gray

Dura-Lite Composite Work Surface Sizes		
Wide	Deep	Cat. No.
36"	30"	20313
48"	30"	20413
60"	30"	20513
72"	30"	20613
96"	30"	20813

Visit [www.HEMCOcorp.com](http://www.HEMCOcorp.com)