


Micro Flow I Ductless Workstation



The MicroFlow I Workstation is ideal for: A Class 1 ductless hood workstation. Activated Carbon Filtration for fumes, odors, and non-hazardous chemical vapors.



- Completely self contained with integral recessed work surface to contain spills. Clear hood surround with safety viewing sash for user. Sash can be conformed for use with a microscope.
- Variable speed fan control. High speed for 100f/m. Air flow thru sash opening. Medium and Low flow for sensitive operations.
- Operates on 115v AC, or 230v International, Conforms to UL, CSA, & CE requirements. Electrical cord port exits left side.

Typical Applications:

- General chemistry involving small volumes of common chemicals, Individual Work stations
- Tissue staining and processing
- Gluing & Drying operations
- Solvent cleaning of electronic parts
- Soldering fumes and odors
- School demonstration workstations
- Containment of Forensic applications



Replacement of pre-filter on Micro Flow I

Specifications:

1. 115v. MicroFlow I Workstation Model # 54010
2. Replacement Carbon Filter Model # 52002
3. Replacement Pre-Filter Model # 54011

Visit www.HEMCOcorp.com