

Multi Purpose Safety Station

The Multi Purpose Safety Station features:



1. Eye/Face Wash
2. Overhead Drench Shower
3. Safety Storage Cabinet

- All in only 4, 5 or 6 sq. ft. of Floor Space!
- Constructed of 3/8" Polypro for durability, sanitation & corrosion resistance.
- Shipped completely assembled, ready to install.
- All plumbing is ready to connect.
- Fixtures meet ANSI Z358.1 standard for Safety Showers.
- Ideal for work sites where persons may be exposed to potentially hazardous chemicals.
- The HEMCO Safety Station is equipped with a pull rod activated shower and push handle eye/face wash for quick rinsing of eyes, face and body.
- Towels and first aid supplies can be stored in the cabinet below.
- Safety Stations are fully assembled and ready for installation to water supply and drain.
- Supply 1.5" NPT, Drain 1.25" NPT.



Safety Station Specifications

"A" Width

A = 24" wide x 84" high x 24" deep. Unit is equipped with Drench Shower, Eye/Face Wash & Storage Cabinet


Product No. 78000

A = 30" wide x 84" high x 24" deep. Unit is equipped with Drench Shower, Eye/Face Wash & Storage Cabinet

Product No. 78001

A = 36" wide x 84" high x 24" deep. Unit is equipped with Drench Shower, Eye/Face Wash & Storage Cabinet

Product No. 78002



Visit www.HEMCOcorp.com

HEMCO Corporation
711 S. Powell Road
Independence, MO 64056

HEMCO

Phone (816) 796-2900
Fax (816) 796-3333
SS-11.11