


Ductless Demonstration Hood

The HEMCO Ductless Demo Hood effectively removes fumes and odors while providing a clear, unobstructed view of the interior from all sides. The demonstration hood is designed for classroom and group observation situations.



- A Class 1 ductless hood workstation.
- Dimensions 35" wide X 24" deep X 30.5 " high
- Viewing height is 24" through clear acrylic, shatterproof side and rear panels.
- Includes hinged viewing sash, phenolic resin work surface, 158 CFM @ 4" speed fan, and power cord.
- Equipped with a vapor proof fluorescent light.
- Front panel includes light switch, variable speed blower, switch and filter service light.
- Activated Carbon Filtration for fumes, odors, non-hazardous chemical vapors.



Specifications:

Demo Hood: Cat. No. 53015 115 V AC, 60 Hz

Mobile Table: Cat. No 50318 37" X 25"
Height adjustable from 30" to 36"

Replacement Filter: Cat. No. 52005 Qty 2 pack

Shipping Weight: 250 lbs. or 113.64 kilos.

Visit www.HEMCOcorp.com

