

# Acrylic Vented Enclosure



## Acrylic Vented Enclosure

### Cat. No. 92200

Ideal size for micro balances. 22" width x 18" depth x 16" height . Enclosure has a 25° ergonomic viewing sash that swings up to a full 15" high access opening. A 6" high reach in height opening is below the sash. Constructed of 3/8 thick clear acrylic with slotted rear baffle and 4" dia. vent outlet.

## Clean-Aire Filtration Module I

### Cat. No. 50110

Convenient solution to point of use exhaust and filtration. Module effectively collects particulates and/or fumes. Can locate on bench top, floor, or shelf and connects to enclosure with flex ducting. Unit operates on 115v and is equipped with variable speed control from 0-250 CFM and 8' power cord. Note: Filters are not included with module.

*Please select filters required below.  
We recommend ordering replacement backup filters.*

Prefilter 9¼" x 8¼" Cat. No. 50111

HEPA Filter 95% , 9¼" x 8¼" Cat. No. 50112

Tri-Media Carbon Filter 9¼" x 8¼" Cat. No. 50113

Acrylic Vented Enclosure - Cat. No. 92200 shown with optional C-Frame Table, Phenolic Worksurface, Flex Ducting and Clean-Aire Filtration Module I.

DESCRIPTION	Cat. No.
<b>Acrylic Vented Enclosure</b> 22" W x 18" D x 16" H with 4" Outlet	92200
<b>Phenolic Worksurface</b> 24" W x 24" D x 1" Thick	12418
<b>C-Frame Table</b> 30" Sitting Height	51203
<b>C-Frame Table</b> 30" Sitting Height w/ casters	51213
<b>Flex Ducting</b> 4" Dia. x 12' long w/ clamps	80051
<b>Exhaust Fan</b> 4" Dia. Inlet /Outlet. 135 CFM	51052
<b>Vent Outlet (thru-wall)</b> 4" Dia. x 15" long	80074



Visit [www.HEMCOcorp.com](http://www.HEMCOcorp.com)

HEMCO Corporation  
711 S. Powell Road  
Independence, MO 64056

# HEMCO

Phone (816) 796-2900  
Fax (816) 796-3333  
AVE-08.10