

# Vented Demonstration Hood

*The HEMCO Ductless Demo Hood effectively removes fumes and odors while providing a clear, unobstructed view of the interior from all sides. The demonstration hood is designed for classroom and group observation situations.*



- Model 53001 Demonstration Hood
- Dimensions 35"wide X 24"deep X 30"high
- Viewing height is 24" through the clear acrylic, shatterproof side and rear panels. Experiments can be safely viewed from the sides and back.
- Sash can be a hinged viewing sash or horizontal sliding sash, both allow easy access.
- The motor housing is chemical resistant, flame retardant fiberglass, the integral motor blower with front mounted switch is designed to be 115 volt/ 60Hertz and provide 12 CFM at 0.4" static pressure.
- Must be connected to an exhaust system.
- Work access area should be 80-100 fpm face velocity.

## Accessories:



### Quick Connect

3inch diameter. Connects exhaust duct to the flexible ductings. Connect is constructed out of polyglass with swing lock arms for quick release and a tight seal **Cat. No. 50130**



### Flexible Ducting

Constructed from PVC material, includes stainless steel clamps. Shipped in 12' lengths. **Cat. No. 80047**



### PVC Vent Outlet

Outlet is for horizontal venting through a wall or window. Angled end prevents rain from entering. male connector is included. Connects to 3" flex ducting. **Cat. No. 80073**

## Optional Table:



### Mobile Table

The table is adjustable from 30 to 36 inches in height. Is of welded steel construction with a plastic laminate work surface. Table comes complete with swivel casters, two of which are lockable. **Cat. No. 50318**



Visit [www.HEMCOcorp.com](http://www.HEMCOcorp.com)